

Tutorial : Control of a C-Model Gripper Using the Modbus RTU Protocol

Description: This tutorial explains how to use the "robotiq_c_model_control" and "robotiq_modbus_rtu" packages to control a C-Model Gripper configured with the Modbus TCP protocol. The C-Model Gripper has 2 fingers. Please visit this website for more information on the [Robotiq Adaptive Robot Grippers](#).

Tutorial level: Beginner

1 - Prerequisites

This tutorial assumes that you have connected a C-Model Gripper via the serial connection with a serial to USB converter to a USB port of your computer. For more information on the Gripper installation please look at the [Robotiq C-Model user manual](#). If you wish to connect the Gripper via the TCP protocol, please have a look at the tutorial on Control of an S-Model Gripper using the TCP Protocol. To test your Gripper and the serial interface you can try the windows driver as explained in the [Robotiq software installation page](#).

Finally, make sure the external dependency for the package "robotiq_modbus_tcp" has been installed. The dependency is the python package pyModbus. On Ubuntu Precise (12.04), it is simply installed using:

```
$ rosdep install robotiq_modbus_tcp
```

On other systems, it can be installed using:

```
$ easy_install -U pymodbus
```

It is possible that you are required to have administrator rights. If so, simply add "sudo" at the beginning of the command.

2 - Configuring the Serial Port

To find out the port on which the controller is connected, use:

```
$ dmesg | grep tty
```

It's normally attributed to ttyUSB0, which we can then call with /dev/ttyUSB0. You might need to change the access permissions on the USB port. It can be done with:

```
$ sudo chmod 777 /dev/ttyUSB0
```

Once this is done, the driver's node can be started.

3 - Starting the Node

3.1 - Run the C-Model Driver Node

After sending the usual "roscore" command, the Gripper can be driven by the node "CModelRtuNode.py" contained in the package "robotiq_c_model_control". The USB port name has to be provided as an argument.

For example, the driver for controlling a C-Model Gripper with the USB port name /dev/ttyUSB0 is launched with the following command:

```
$ rosrun robotiq_c_model_control CModelRtuNode.py /dev/ttyUSB0
```

When the Gripper is successfully connected, the LED light appears blue.

3.2 - Run the C-Model Simple Controller Node

The driver listens for messages on "CModelRobotOutput" using the "SModel_robot_output" msg type. The messages are interpreted and commands are sent to the Gripper accordingly. A simple controller node is provided which can be run in another terminal, using the command:

```
$ rosrun robotiq_c_model_control CModelSimpleController.py
```

```
-----  
Current command:  rACT = 0, rGTO = 0, rATR = 0, rPR = 0, rSP = 0, rFR = 0  
-----  
Available commands  
  
r: Reset  
a: Activate  
c: Close  
o: Open  
(0-255): Go to that position  
f: Faster  
l: Slower  
i: Increase force  
d: Decrease force  
-->
```

Figure 1 - Simple Controller Node

The "CModel_robot_output" msg type is simply composed of the robot output variables described in the [Robotiq C-Model user manual](#). The simple controller node can therefore be modified to send custom commands to the Gripper.

3.3 - Run the C-Model Status Listener Node

In the package "robotiq_c_model_control", there is also a node for listening and interpreting the status of the Gripper. The driver publishes the status of the Gripper on "CModelRobotInput" using the "CModel_robot_input" msg type. The msg type is composed of the robot input variables described in the [Robotiq C-Model user manual](#). The status listener node can be run in another terminal, using the following command:

```
$ rosrun robotiq_c_model_control CModelStatusListener.py
```

```
-----  
C-Model status interpreter  
-----  
gACT = 1: Gripper activation  
gGTO = 1: Go to Position Request  
gSTA = 3: Activation is completed  
gOBJ = 3: Fingers are at requested position  
gFLT = 0: No Fault  
gPR = 0: Echo of the requested position for the Gripper: 0/255  
gPO = 3: Position of Fingers: 3/255  
gCU = 0: Current of Fingers: 0 mA
```

Figure 2 - C-Model Status Listener